

Integral Version 3.0 Announcement -20180806

Labsphere, the world leader in light measurement, is proud to announce our latest **illumia Plus** software release: **Integral 3.0**

Integral 3.0 adds significant new features to **illumia Plus** since V2.0. These changes improve the quality, user experience and functionality of the software and provide greater flexibility, added analysis capability and user-friendly features, all designed to enhance your overall experience.

NEW FEATURES IN INTEGRAL 3.0

Temporal Lighting Artifacts¹! This new feature with the Flicker-iP accessory include

- Selectable Sampling Rate
- Selectable Scan Duration (Measurement Period)
- Selectable Recording Intervals
- Lamp Light Output Periodic Frequency
- Percent Flicker
- Flicker Index
- Stroboscopic Visibility Measure (SVM)
- Short Term Flicker (P_{stLM})
- ASSIST Flicker Perception Metric (M_p)
- Fraction of Rated Light Output Integrated over Measurement Period
- Percent Amplitude Modulation at Selected Frequencies
- Precision Cosine Receiver
- Robust Carry and Storage Case

Project Scan List Enhancements

- Search Bar and Reset
- Customization of columns
- Multi-select and deselect checkboxes from anchor point
- Scanlist height animation when new scans are saved

YOUR UPGRADE INFORMATION

Labsphere's goal is to make the software upgrade process as seamless and effortless as possible for you, our valued customers.

Just go to the **direct link to Integral releases and updates from your Integral Menu**, or access registration for V3.0 by following the link below

<https://www.labsphere.com/support/software-downloads-private/>

Complete the registration form and a Customer Service Representative will contact you with further instructions.

¹ Available with Integral 3.0 or Purchase of Flicker-iP

Integral Version 3.0 Announcement -20180806

We look forward to our continued partnership with you and your company. If you have any questions, concerns or need help, please request Labsphere's Service Support by +1(603) 927-4266 or by e-mail at Labsphere@labsphere.com

We appreciate any feedback you can offer us on **Integral** software or any other Labsphere products or services.